

Product datasheet for **SR421937**

Sall4 Mouse siRNA Oligo Duplex (Locus ID 99377)

Product data:

Product Type:	siRNA Oligo Duplexes
Purity:	HPLC purified
Quality Control:	Tested by ESI-MS
Sequences:	Available with shipment
Stability:	One year from date of shipment when stored at -20°C.
# of transfections:	Approximately 330 transfections/2nmol in 24-well plate under optimized conditions (final conc. 10 nM).
Note:	Single siRNA duplex (10nmol) can be ordered.
RefSeq:	NM_175303 , NM_201395 , NM_201396
UniProt ID:	Q8BX22
Synonyms:	5730441M18Rik; AA407717; AL022809; AW536104; C78083; C78563; C330011P20Rik; Tex20
Components:	Sall4 (Mouse) - 3 unique 27mer siRNA duplexes - 2 nmol each (Locus ID 99377) Included - SR30004, Trilencer-27 Universal Scrambled Negative Control siRNA Duplex - 2 nmol Included - SR30005, RNase free siRNA Duplex Resuspension Buffer - 2 ml
Summary:	This gene belongs to the spalt family of zinc finger transcription factors. In mouse, functions for this gene have been described in many embryonic developmental processes, including brain, heart, and limb development. In addition, this gene is an important pluripotency factor that is required for stem cell maintenance. Homozygous mutant mice display embryonic lethality, while conditional knock-out in embryonic germ cells results in failure to establish a robust stem cell population. A pseudogene of this gene is found on chromosome 2. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Apr 2015]



[View online »](#)

**Performance
Guaranteed:**

OriGene guarantees that at least two of the three Dicer-Substrate duplexes in the kit will provide at least 70% or more knockdown of the target mRNA when used at 10 nM concentration by quantitative RT-PCR when the TYE-563 fluorescent transfection control duplex (cat# SR30002) indicates that >90% of the cells have been transfected and the HPRT positive control (cat# SR30003) provides 90% knockdown efficiency.

For non-conforming siRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the siRNA kit. To arrange for a free replacement with newly designed duplexes, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled siRNA control (quantitative RT-PCR data required).