

## Product datasheet for **SR422158**

### Atp11c Mouse siRNA Oligo Duplex (Locus ID 320940)

#### Product data:

Product Type:	siRNA Oligo Duplexes
Purity:	HPLC purified
Quality Control:	Tested by ESI-MS
Sequences:	Available with shipment
Stability:	One year from date of shipment when stored at -20°C.
# of transfections:	Approximately 330 transfections/2nmol in 24-well plate under optimized conditions (final conc. 10 nM).
Note:	Single siRNA duplex (10nmol) can be ordered.
RefSeq:	<a href="#">NM_001001798</a> , <a href="#">NM_001037863</a> , <a href="#">NM_001359002</a>
UniProt ID:	<a href="#">Q9QZW0</a>
Synonyms:	A330005H02Rik; AI315324; Ig
Components:	Atp11c (Mouse) - 3 unique 27mer siRNA duplexes - 2 nmol each (Locus ID 320940) Included - SR30004, Trilencer-27 Universal Scrambled Negative Control siRNA Duplex - 2 nmol Included - SR30005, RNase free siRNA Duplex Resuspension Buffer - 2 ml
Summary:	Catalytic component of a P4-ATPase flippase complex which catalyzes the hydrolysis of ATP coupled to the transport of aminophospholipids from the outer to the inner leaflet of various membranes and ensures the maintenance of asymmetric distribution of phospholipids. In the cell membrane of erythrocytes, it is required to maintain phosphatidylserine (PS) in the inner leaflet preventing its exposure on the surface. This asymmetric distribution is critical for the survival of erythrocytes in circulation since externalized PS is a phagocytic signal for splenic macrophages (By similarity). Phospholipid translocation seems also to be implicated in vesicle formation and in uptake of lipid signaling molecules. Required for B cell differentiation past the pro-B cell stage (PubMed:21423173). Seems to mediate phosphatidylserine (PS) flipping in pro-B cells (PubMed:21423172). May be involved in the transport of cholestatic bile acids (PubMed:21518881).[UniProtKB/Swiss-Prot Function]


[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that at least two of the three Dicer-Substrate duplexes in the kit will provide at least 70% or more knockdown of the target mRNA when used at 10 nM concentration by quantitative RT-PCR when the TYE-563 fluorescent transfection control duplex (cat# SR30002) indicates that >90% of the cells have been transfected and the HPRT positive control (cat# SR30003) provides 90% knockdown efficiency.

For non-conforming siRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the siRNA kit. To arrange for a free replacement with newly designed duplexes, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled siRNA control (quantitative RT-PCR data required).