

Product datasheet for **SR422505**

Ctnnd2 Mouse siRNA Oligo Duplex (Locus ID 18163)

Product data:

Product Type:	siRNA Oligo Duplexes
Purity:	HPLC purified
Quality Control:	Tested by ESI-MS
Sequences:	Available with shipment
Stability:	One year from date of shipment when stored at -20°C.
# of transfections:	Approximately 330 transfections/2nmol in 24-well plate under optimized conditions (final conc. 10 nM).
Note:	Single siRNA duplex (10nmol) can be ordered.
RefSeq:	NM_008729 , NM_001355502
UniProt ID:	Q35927
Synonyms:	Catnd2; neurojugin; Nprap
Components:	Ctnnd2 (Mouse) - 3 unique 27mer siRNA duplexes - 2 nmol each (Locus ID 18163) Included - SR30004, Trilencer-27 Universal Scrambled Negative Control siRNA Duplex - 2 nmol Included - SR30005, RNase free siRNA Duplex Resuspension Buffer - 2 ml
Summary:	Has a critical role in neuronal development, particularly in the formation and/or maintenance of dendritic spines and synapses (PubMed:17993462) (PubMed:25807484). Involved in the regulation of canonical Wnt signaling (By similarity). It probably acts on beta-catenin turnover, facilitating beta-catenin interaction with GSK3B, phosphorylation, ubiquitination and degradation (PubMed:20623542). May be involved in neuronal cell adhesion and tissue morphogenesis and integrity by regulating adhesion molecules. Functions as a transcriptional activator when bound to ZBTB33 (PubMed:15282317).[UniProtKB/Swiss-Prot Function]



[View online »](#)

**Performance
Guaranteed:**

OriGene guarantees that at least two of the three Dicer-Substrate duplexes in the kit will provide at least 70% or more knockdown of the target mRNA when used at 10 nM concentration by quantitative RT-PCR when the TYE-563 fluorescent transfection control duplex (cat# SR30002) indicates that >90% of the cells have been transfected and the HPRT positive control (cat# SR30003) provides 90% knockdown efficiency.

For non-conforming siRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the siRNA kit. To arrange for a free replacement with newly designed duplexes, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled siRNA control (quantitative RT-PCR data required).