

Product datasheet for **SR422684**

Pard3 Mouse siRNA Oligo Duplex (Locus ID 93742)

Product data:

Product Type:	siRNA Oligo Duplexes
Purity:	HPLC purified
Quality Control:	Tested by ESI-MS
Sequences:	Available with shipment
Stability:	One year from date of shipment when stored at -20°C.
# of transfections:	Approximately 330 transfections/2nmol in 24-well plate under optimized conditions (final conc. 10 nM).
Note:	Single siRNA duplex (10nmol) can be ordered.
RefSeq:	NM_001013580 , NM_001013581 , NM_001122850 , NM_001309391 , NM_001309392 , NM_033620 , NM_001363431 , NM_001363432
UniProt ID:	Q99NH2
Synonyms:	AA960621; AI256638; Asip; D8Ertd580e; Par-3; Par3; Pard-3; Pard3a; Phip
Components:	Pard3 (Mouse) - 3 unique 27mer siRNA duplexes - 2 nmol each (Locus ID 93742) Included - SR30004, Trilencer-27 Universal Scrambled Negative Control siRNA Duplex - 2 nmol Included - SR30005, RNase free siRNA Duplex Resuspension Buffer - 2 ml
Summary:	Adapter protein involved in asymmetrical cell division and cell polarization processes (By similarity). Seems to play a central role in the formation of epithelial tight junctions (By similarity). Targets the phosphatase PTEN to cell junctions (By similarity). Association with PARD6B may prevent the interaction of PARD3 with F11R/JAM1, thereby preventing tight junction assembly (PubMed:11839275). The PARD6-PARD3 complex links GTP-bound Rho small GTPases to atypical protein kinase C proteins (By similarity). Required for establishment of neuronal polarity and normal axon formation in cultured hippocampal neurons (By similarity). Involved in Schwann cell peripheral myelination (PubMed:21949390). [UniProtKB/Swiss-Prot Function]



[View online »](#)

**Performance
Guaranteed:**

OriGene guarantees that at least two of the three Dicer-Substrate duplexes in the kit will provide at least 70% or more knockdown of the target mRNA when used at 10 nM concentration by quantitative RT-PCR when the TYE-563 fluorescent transfection control duplex (cat# SR30002) indicates that >90% of the cells have been transfected and the HPRT positive control (cat# SR30003) provides 90% knockdown efficiency.

For non-conforming siRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the siRNA kit. To arrange for a free replacement with newly designed duplexes, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled siRNA control (quantitative RT-PCR data required).