

## Product datasheet for **SR423005**

### Madd Mouse siRNA Oligo Duplex (Locus ID 228355)

#### Product data:

Product Type:	siRNA Oligo Duplexes
Purity:	HPLC purified
Quality Control:	Tested by ESI-MS
Sequences:	Available with shipment
Stability:	One year from date of shipment when stored at -20°C.
# of transfections:	Approximately 330 transfections/2nmol in 24-well plate under optimized conditions (final conc. 10 nM).
Note:	Single siRNA duplex (10nmol) can be ordered.
RefSeq:	<a href="#">NM_001177719</a> , <a href="#">NM_001177720</a> , <a href="#">NM_001177721</a> , <a href="#">NM_001177722</a> , <a href="#">NM_001177723</a> , <a href="#">NM_001177724</a> , <a href="#">NM_001177725</a> , <a href="#">NM_001177726</a> , <a href="#">NM_001177727</a> , <a href="#">NM_001177728</a> , <a href="#">NM_001177729</a> , <a href="#">NM_145527</a>
UniProt ID:	<a href="#">Q80U28</a>
Synonyms:	9630059K23Rik; IG20
Components:	Madd (Mouse) - 3 unique 27mer siRNA duplexes - 2 nmol each (Locus ID 228355) Included - SR30004, Trilencer-27 Universal Scrambled Negative Control siRNA Duplex - 2 nmol Included - SR30005, RNAse free siRNA Duplex Resuspension Buffer - 2 ml
Summary:	Plays a significant role in regulating cell proliferation, survival and death through alternative mRNA splicing. Converts GDP-bound inactive form of RAB3A, RAB3C and RAB3D to the GTP-bound active forms. Component of the TNFRSF1A signaling complex: MADD links TNFRSF1A with MAP kinase activation. Plays an important regulatory role in physiological cell death (TNF-alpha-induced, caspase-mediated apoptosis).[UniProtKB/Swiss-Prot Function]



[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that at least two of the three Dicer-Substrate duplexes in the kit will provide at least 70% or more knockdown of the target mRNA when used at 10 nM concentration by quantitative RT-PCR when the TYE-563 fluorescent transfection control duplex (cat# SR30002) indicates that >90% of the cells have been transfected and the HPRT positive control (cat# SR30003) provides 90% knockdown efficiency.

For non-conforming siRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the siRNA kit. To arrange for a free replacement with newly designed duplexes, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled siRNA control (quantitative RT-PCR data required).