

## Product datasheet for **SR423215**

### Col11a1 Mouse siRNA Oligo Duplex (Locus ID 12814)

#### Product data:

Product Type:	siRNA Oligo Duplexes
Purity:	HPLC purified
Quality Control:	Tested by ESI-MS
Sequences:	Available with shipment
Stability:	One year from date of shipment when stored at -20°C.
# of transfections:	Approximately 330 transfections/2nmol in 24-well plate under optimized conditions (final conc. 10 nM).
Note:	Single siRNA duplex (10nmol) can be ordered.
RefSeq:	<a href="#">NM_007729</a>
UniProt ID:	<a href="#">Q61245</a>
Synonyms:	C530001D20Rik; cho
Components:	Col11a1 (Mouse) - 3 unique 27mer siRNA duplexes - 2 nmol each (Locus ID 12814) Included - SR30004, Trilencer-27 Universal Scrambled Negative Control siRNA Duplex - 2 nmol Included - SR30005, RNase free siRNA Duplex Resuspension Buffer - 2 ml
Summary:	This gene encodes the alpha-1 subunit of type XI collagen, one of the low abundance fibrillar collagens that is essential for normal embryonic skeletal development and the cohesive properties of cartilage. The encoded protein, in association with the alpha-1 subunit of type II collagen, forms a heterotrimeric type XI procollagen that undergoes proteolytic processing. Mice lacking the encoded protein develop severe chondrodysplasia and die at birth. [provided by RefSeq, Dec 2015]



[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that at least two of the three Dicer-Substrate duplexes in the kit will provide at least 70% or more knockdown of the target mRNA when used at 10 nM concentration by quantitative RT-PCR when the TYE-563 fluorescent transfection control duplex (cat# SR30002) indicates that >90% of the cells have been transfected and the HPRT positive control (cat# SR30003) provides 90% knockdown efficiency.

For non-conforming siRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the siRNA kit. To arrange for a free replacement with newly designed duplexes, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled siRNA control (quantitative RT-PCR data required).