

Product datasheet for SR423518

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Lrrk2 Mouse siRNA Oligo Duplex (Locus ID 66725)

Product data:

UniProt ID:

Product Type: siRNA Oligo Duplexes

Purity: HPLC purified

Quality Control: Tested by ESI-MS

Sequences: Available with shipment

Stability: One year from date of shipment when stored at -20°C.

of transfections: Approximately 330 transfections/2nmol in 24-well plate under optimized conditions (final

conc. 10 nM).

Q5S006

Note: Single siRNA duplex (10nmol) can be ordered.

RefSeq: <u>NM 025730</u>

Synonyms: 4921513O20Rik; 9330188B09Rik; AW561911; cl-46; D630001M17Rik; Gm927

Components: Lrrk2 (Mouse) - 3 unique 27mer siRNA duplexes - 2 nmol each (Locus ID 66725)

Included - SR30004, Trilencer-27 Universal Scrambled Negative Control siRNA Duplex - 2 nmol

Included - SR30005, RNAse free siRNA Duplex Resuspension Buffer - 2 ml

Summary: Positively regulates autophagy through a calcium-dependent activation of the CaMKK/AMPK

signaling pathway. The process involves activation of nicotinic acid adenine dinucleotide phosphate (NAADP) receptors, increase in lysosomal pH, and calcium release from lysosomes. Together with RAB29, plays a role in the retrograde trafficking pathway for recycling proteins, such as mannose 6 phosphate receptor (M6PR), between lysosomes and

the Golgi apparatus in a retromer-dependent manner. Regulates neuronal process

morphology in the intact central nervous system (CNS). Phosphorylates PRDX3. Has GTPase activity (By similarity). Plays an important role in recruiting SEC16A to endoplasmic reticulum

exit sites (ERES) and in regulating ER to Golgi vesicle-mediated transport and ERES

organization (PubMed:25201882).[UniProtKB/Swiss-Prot Function]







Performance Guaranteed:

OriGene guarantees that at least two of the three Dicer-Substrate duplexes in the kit will provide at least 70% or more knockdown of the target mRNA when used at 10 nM concentration by quantitative RT-PCR when the TYE-563 fluorescent transfection control duplex (cat# SR30002) indicates that >90% of the cells have been transfected and the HPRT positive control (cat# SR30003) provides 90% knockdown efficiency.

For non-conforming siRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the siRNA kit. To arrange for a free replacement with newly designed duplexes, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled siRNA control (quantitative RT-PCR data required).