

Product datasheet for **SR426753**

Adrb2 Mouse siRNA Oligo Duplex (Locus ID 11555)

Product data:

Product Type:	siRNA Oligo Duplexes
Purity:	HPLC purified
Quality Control:	Tested by ESI-MS
Sequences:	Available with shipment
Stability:	One year from date of shipment when stored at -20°C.
# of transfections:	Approximately 330 transfections/2nmol in 24-well plate under optimized conditions (final conc. 10 nM).
Note:	Single siRNA duplex (10nmol) can be ordered.
RefSeq:	NM_007420
UniProt ID:	P18762
Synonyms:	Adrb; Adrb-2; Bad; Badm; Gpc; Gpcr7
Components:	Adrb2 (Mouse) - 3 unique 27mer siRNA duplexes - 2 nmol each (Locus ID 11555) Included - SR30004, Trilencer-27 Universal Scrambled Negative Control siRNA Duplex - 2 nmol Included - SR30005, RNase free siRNA Duplex Resuspension Buffer - 2 ml
Summary:	<p>This intronless gene belongs to the G-protein-coupled receptor superfamily, which includes transmembrane proteins that play a role in signal transduction across biological membranes resulting in a variety of physiological responses. The encoded protein is a beta-2 adrenergic receptor which is activated by catecholamine ligands such as adrenaline and epinephrine. The protein participates in the classical signaling pathway involving G protein, adenylyl cyclase, cAMP (3'-5'-cyclic adenosine monophosphate) and protein kinase A (PKA). In humans, this gene is implicated in susceptibility to nocturnal asthma, obesity and type 2 diabetes. [provided by RefSeq, Apr 2013]</p>



[View online »](#)

**Performance
Guaranteed:**

OriGene guarantees that at least two of the three Dicer-Substrate duplexes in the kit will provide at least 70% or more knockdown of the target mRNA when used at 10 nM concentration by quantitative RT-PCR when the TYE-563 fluorescent transfection control duplex (cat# SR30002) indicates that >90% of the cells have been transfected and the HPRT positive control (cat# SR30003) provides 90% knockdown efficiency.

For non-conforming siRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the siRNA kit. To arrange for a free replacement with newly designed duplexes, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled siRNA control (quantitative RT-PCR data required).