

Product datasheet for SR506650

Park7 Rat siRNA Oligo Duplex (Locus ID 117287)

Product data:

Product Type: siRNA Oligo Duplexes HPLC purified **Purity: Quality Control:** Tested by ESI-MS Available with shipment Sequences: Stability: One year from date of shipment when stored at -20°C. # of transfections: Approximately 330 transfections/2nmol in 24-well plate under optimized conditions (final conc. 10 nM). Note: Single siRNA duplex (10nmol) can be ordered. **RefSeq:** NM 001277249, NM 001277250, NM 001277251, NM 001277252, NM 001277253, NM 057143 UniProt ID: 088767 Synonyms: CAP1; DJ-1; Dj1; SP22 Park7 (Rat) - 3 unique 27mer siRNA duplexes - 2 nmol each (Locus ID 117287) **Components:** Included - SR30004, Trilencer-27 Universal Scrambled Negative Control siRNA Duplex - 2 nmol Included - SR30005, RNAse free siRNA Duplex Resuspension Buffer - 2 ml displays a specific localization pattern to the anterior-ventral surface of the equatorial Summary: segment; may play a role in sperm-egg interactions and fertilization [RGD, Feb 2006] Performance OriGene guarantees that at least two of the three Dicer-Substrate duplexes in the kit will **Guaranteed:** provide at least 70% or more knockdown of the target mRNA when used at 10 nM concentration by quantitative RT-PCR when the TYE-563 fluorescent transfection control duplex (cat# SR30002) indicates that >90% of the cells have been transfected and the HPRT positive control (cat# SR30003) provides 90% knockdown efficiency. For non-conforming siRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the siRNA kit. To arrange for a free replacement with newly designed duplexes, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene

expression knockdown compared to the scrambled siRNA control (quantitative RT-PCR data



/iew online »

This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2021 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US

required).

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn