

Product datasheet for **SR508475**

Nrxn1 Rat siRNA Oligo Duplex (Locus ID 60391)

Product data:

Product Type:	siRNA Oligo Duplexes
Purity:	HPLC purified
Quality Control:	Tested by ESI-MS
Sequences:	Available with shipment
Stability:	One year from date of shipment when stored at -20°C.
# of transfections:	Approximately 330 transfections/2nmol in 24-well plate under optimized conditions (final conc. 10 nM).
Note:	Single siRNA duplex (10nmol) can be ordered.
RefSeq:	NM_021767
UniProt ID:	Q63372
Synonyms:	Nrxn1b
Components:	Nrxn1 (Rat) - 3 unique 27mer siRNA duplexes - 2 nmol each (Locus ID 60391) Included - SR30004, Trilencer-27 Universal Scrambled Negative Control siRNA Duplex - 2 nmol Included - SR30005, RNase free siRNA Duplex Resuspension Buffer - 2 ml
Summary:	This gene encodes a single-pass type I membrane protein that belongs to the neurexin family. Neurexins are synaptic transmembrane receptors that bind endogenous ligands that include neuroligins, dystroglycan, and neurexophilins. Neurexin complexes are required for efficient neurotransmission and are involved in synaptogenesis. In vertebrates, alternate promoter usage results in multiple isoform classes, of which the alpha and beta classes are the best characterized. In humans, allelic variants in this gene are associated with Pitt-Hopkins-like syndrome-2, while deletions have been associated with autism and schizophrenia. Mouse knockouts display decreased spontaneous and evoked vesicle release resulting in impaired synaptic transmission. In addition, knockout mice show altered social approach, reduced social investigation, reduced locomotor activity, and in males, increased aggression. [provided by RefSeq, Nov 2016]



[View online »](#)

**Performance
Guaranteed:**

OriGene guarantees that at least two of the three Dicer-Substrate duplexes in the kit will provide at least 70% or more knockdown of the target mRNA when used at 10 nM concentration by quantitative RT-PCR when the TYE-563 fluorescent transfection control duplex (cat# SR30002) indicates that >90% of the cells have been transfected and the HPRT positive control (cat# SR30003) provides 90% knockdown efficiency.

For non-conforming siRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the siRNA kit. To arrange for a free replacement with newly designed duplexes, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled siRNA control (quantitative RT-PCR data required).