

## Product datasheet for **SR514355**

### Gskip Rat siRNA Oligo Duplex (Locus ID 362778)

#### Product data:

Product Type:	siRNA Oligo Duplexes
Purity:	HPLC purified
Quality Control:	Tested by ESI-MS
Sequences:	Available with shipment
Stability:	One year from date of shipment when stored at -20°C.
# of transfections:	Approximately 330 transfections/2nmol in 24-well plate under optimized conditions (final conc. 10 nM).
Note:	Single siRNA duplex (10nmol) can be ordered.
RefSeq:	<a href="#">NM_001014198</a>
UniProt ID:	<a href="#">Q5PPI3</a>
Synonyms:	RGD1308470
Components:	Gskip (Rat) - 3 unique 27mer siRNA duplexes - 2 nmol each (Locus ID 362778) Included - SR30004, Trilencer-27 Universal Scrambled Negative Control siRNA Duplex - 2 nmol Included - SR30005, RNase free siRNA Duplex Resuspension Buffer - 2 ml
Summary:	A-kinase anchoring protein for GSK3B and PKA that regulates or facilitates their kinase activity towards their targets. The ternary complex enhances Wnt-induced signaling by facilitating the GSK3B- and PKA-induced phosphorylation of beta-catenin leading to beta-catenin degradation and stabilization respectively. Upon cAMP activation, the ternary complex contributes to neuroprotection against oxidative stress-induced apoptosis by facilitating the PKA-induced phosphorylation of DML1 and PKA-induced inactivation of GSK3B. During neurite outgrowth promotes neuron proliferation; while increases beta-catenin-induced transcriptional activity through GSK3B kinase activity inhibition, reduces N-cadherin level to promote cell cycle progression (By similarity). May play a role in cleft palate formation and is required for postnatal life through modulation of the activity of GSK3B during development (By similarity).[UniProtKB/Swiss-Prot Function]



[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that at least two of the three Dicer-Substrate duplexes in the kit will provide at least 70% or more knockdown of the target mRNA when used at 10 nM concentration by quantitative RT-PCR when the TYE-563 fluorescent transfection control duplex (cat# SR30002) indicates that >90% of the cells have been transfected and the HPRT positive control (cat# SR30003) provides 90% knockdown efficiency.

For non-conforming siRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the siRNA kit. To arrange for a free replacement with newly designed duplexes, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled siRNA control (quantitative RT-PCR data required).