

## Product datasheet for **TA100031**

### Rabbit Polyclonal Anti-GST Tag Antibody


#### Product data:

Product Type:	Tag Antibodies
Product Name:	Rabbit Polyclonal Anti-GST Tag Antibody
Applications:	WB
Recommended Dilution:	WB: 2000-10000 WB positive control: Recombinant GST tag protein
Reactivity:	Fusion proteins containing the GST Tag
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	Full length fusion protein
Isotype:	IgG
Formulation:	pH7.4 PBS, 0.05% NaN <sub>3</sub> , 40% Glycerol
Concentration:	2.5 mg/ml
Conjugation:	Unconjugated
Storage:	Store at -20°C. For maximum recovery of product, centrifuge the original vial after thawing and prior to opening the cap. Aliquot to avoid repeated freezing and thawing.
Stability:	Stable for 1 year at -20°C from date of shipment.
Predicted Protein Size:	63 kDa



[View online »](#)

## Product images:



Gel: 8%SDS-PAGE  
Lysate: 0.005/0.01/0.05/0.1  $\mu$ g  
Lane 1-4: 0.005/0.01/0.05/0.1  $\mu$ g Recombinant  
GST tag protein  
Primary antibody: TA100031 (GST Tag Antibody)  
at dilution 1/4000  
Secondary antibody: Goat anti rabbit IgG at  
1/5000 dilution  
Exposure time: 30 seconds