

Product datasheet for **TA130002**

Rabbit Anti-Mouse IgG Secondary Antibody

Product data:

Product Type:	Secondary Antibodies
Product Name:	Rabbit Anti-Mouse IgG Secondary Antibody
Applications:	ELISA, IHC, WB
Recommended Dilution:	WB, IHC, IF/ICC, ELISA
Reactivity:	Mouse IgG
Host:	Rabbit
Immunogen:	Mouse IgG, whole molecule
Isotype:	IgG
Formulation:	Lyophilized with 1% sodium acetate
Purification:	This antibody was purified from antiserum by immunoaffinity chromatography which removes essentially all rabbit serum proteins, except the specific antibody for mouse IgG.
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.



[View online »](#)