

Product datasheet for **TA130004**

Peroxidase Conjugated Goat Anti-Mouse IgG (gamma-chain specific)

Product data:

Product Type:	Secondary Antibodies
Product Name:	Peroxidase Conjugated Goat Anti-Mouse IgG (gamma-chain specific)
Applications:	Dot
Recommended Dilution:	Dot blot:0.5-1µg/ml.Western blot(DAB):0.7-3.3µg/ml.Western blot(ECL):0.2-0.5µg/ml.ELISA:0.2-0.3µg/ml.
Reactivity:	Mouse IgG
Host:	Goat
Immunogen:	Mouse IgG, whole molecule
Isotype:	IgG
Formulation:	0.01M PBS (pH7.4), 50% glycerol
Concentration:	lot specific
Purification:	This antibody was purified from antiserum by immunoaffinity chromatography which removes essentially all goat serum proteins, except the specific antibody for mouse IgG.
Conjugation:	HRP
Storage:	Store at +4°C for one year.
Note:	This peroxidase conjugated antibody is specific for mouse IgG and shows no cross-reactivity with mouse IgA, IgM and other animal species IgG proteins.



[View online »](#)