

## Product datasheet for **TA130005**

### Peroxidase Conjugated Goat Anti-Mouse IgG (gamma-chain specific)

#### Product data:

Product Type:	Secondary Antibodies
Product Name:	Peroxidase Conjugated Goat Anti-Mouse IgG (gamma-chain specific)
Applications:	ELISA, WB
Recommended Dilution:	WB: 1:3,000-10,000; ELISA: 1:20,000-40,000
Reactivity:	Mouse IgG
Host:	Goat
Immunogen:	Mouse IgG, whole molecule
Isotype:	IgG
Formulation:	0.01M PBS (pH7.4), 50% glycerol
Concentration:	lot specific
Purification:	This antibody was purified from antiserum by immunoaffinity chromatography which removes essentially all goat serum proteins, except the specific antibody for mouse IgG.
Conjugation:	HRP
Storage:	Store at -20°C as received. Avoid repeated freezing and thawing. After first thaw, pipette into single aliquote for future use.
Note:	This peroxidase conjugated antibody is specific for mouse IgG and shows no cross-reactivity with mouse IgA, IgM and other animal species IgG proteins.



[View online »](#)