

Product datasheet for TA130009

Biotin Conjugated Goat Anti-Mouse IgG

Product data:

Product Type:	Secondary Antibodies
Product Name:	Biotin Conjugated Goat Anti-Mouse IgG
Applications:	ELISA, ICC, IHC, WB
Recommend Dilution:	WB: 1:200-400; ICC: 1:100-200; IHC: 1:100-200; ELISA: 1:5,000-20,000
Reactivity:	Mouse IgG
Host:	Goat
Immunogen:	Mouse IgG, whole molecule
Isotype:	IgG
Formulation:	0.01M PBS (pH7.4), 0.01% Thimerosal, 50% glycerol
Form:	Concentrated, Liquid
Concentration:	lot specific
Purification:	This antibody was purified from antiserum by immunoaffinity chromatography which removes essentially all goat serum proteins, except the specific antibody for mouse IgG.
Conjugation:	Biotin
Storage:	Store at -20°C as received. Avoid repeated freezing and thawing. After first thaw, pipette into single aliquote for future use.
Note:	This biotin conjugated antibody is specific for mouse IgG and shows no cross-reactivity with mouse IgA, IgM and other animal species IgG proteins.



[View online »](#)