

## Product datasheet for **TA130010**

### Biotin Conjugated Goat Anti-Mouse IgM

#### Product data:

Product Type:	Secondary Antibodies
Product Name:	Biotin Conjugated Goat Anti-Mouse IgM
Applications:	ELISA, ICC, IHC, WB
Recommended Dilution:	WB: 1:200-400; ICC: 1:100-200; IHC: 1:100-200; ELISA: 1:5,000-20,000
Reactivity:	Mouse IgM
Host:	Goat
Immunogen:	Mouse IgM, whole molecule
Isotype:	IgM
Formulation:	0.01M PBS (pH7.4), 0.01% Thimerosal, 50% glycerol
Concentration:	lot specific
Purification:	This antibody was purified from antiserum by immunoaffinity chromatography which removes essentially all goat serum proteins, except the specific antibody for mouse IgM.
Conjugation:	Biotin
Storage:	Store at -20°C as received. Avoid repeated freezing and thawing. After first thaw, pipette into single aliquote for future use.
Note:	This biotin conjugated antibody is specific for mouse IgM and shows no cross-reactivity with mouse IgA, IgG and other animal species IgM proteins.



[View online »](#)