

## Product datasheet for TA130012

### Cy3 Conjugated Goat Anti-Mouse IgG

#### Product data:

Product Type:	Secondary Antibodies
Product Name:	Cy3 Conjugated Goat Anti-Mouse IgG
Applications:	IHC
Recommend Dilution:	IF/ICC: 1:50-100; IHC: 1:50-100
Reactivity:	Mouse IgG
Host:	Goat
Immunogen:	Mouse IgG, whole molecule
Isotype:	IgG
Formulation:	0.01M PBS (pH7.4), 0.01% Thimerosal, 50% glycerol
Form:	Concentrated, Liquid
Concentration:	lot specific
Purification:	This antibody was purified from antiserum by immunoaffinity chromatography which removes essentially all goat serum proteins, except the specific antibody for mouse IgG.
Conjugation:	Cy3
Storage:	Store at -20°C as received. Avoid repeated freezing and thawing. After first thaw, pipette into single aliquote for future use.



[View online »](#)