

## Product datasheet for **TA130014**

### **FITC Conjugated Goat Anti-Mouse IgG**

#### **Product data:**

|                       |   |
|-----------------------|---|
| Product Type:         | Secondary Antibodies  |
| Product Name:         | FITC Conjugated Goat Anti-Mouse IgG   |
| Applications:         | FC, IHC   |
| Recommended Dilution: | FC: 1:100; IF/ICC: 1:32-64; IHC: 1:32-64  |
| Reactivity:           | Mouse IgG   |
| Host:                 | Goat  |
| Immunogen:            | Mouse IgG, whole molecule   |
| Isotype:              | IgG   |
| Formulation:          | 0.01M PBS (pH7.4), 0.01% Thimerosal, 50% glycerol   |
| Concentration:        | lot specific  |
| Purification:         | This antibody was purified from antiserum by immunoaffinity chromatography which removes essentially all goat serum proteins, except the specific antibody for mouse IgG. |
| Conjugation:          | FITC  |
| Storage:              | Store at -20°C as received. Avoid repeated freezing and thawing. After first thaw, pipette into single aliquote for future use.   |
| Note:                 | This FITC conjugated antibody is specific for mouse IgG and shows no cross-reactivity with mouse IgA, IgM and other animal species IgG proteins.                          |



[View online »](#)