

## Product datasheet for **TA130015**

### Goat Anti-Rabbit IgG Secondary Antibody

#### Product data:

Product Type:	Secondary Antibodies
Product Name:	Goat Anti-Rabbit IgG Secondary Antibody
Applications:	ELISA, IHC, WB
Recommended Dilution:	WB, IHC, IF/ICC, ELISA
Reactivity:	Rabbit IgG
Host:	Goat
Immunogen:	Rabbit IgG, whole molecule
Isotype:	IgG
Formulation:	Lyophilized with 1% sodium acetate
Purification:	This antibody was purified from antiserum by immunoaffinity chromatography which removes essentially all goat serum proteins, except the specific antibody for rabbit IgG.
Conjugation:	Unconjugated
Storage:	Store at -20°C as received. Avoid repeated freezing and thawing. After first thaw, pipette into single aliquote for future use.



[View online »](#)