

Product datasheet for **TA130019**

Biotin Conjugated Mouse Anti-Rabbit IgG

Product data:

Product Type:	Secondary Antibodies
Product Name:	Biotin Conjugated Mouse Anti-Rabbit IgG
Applications:	ELISA, ICC, IHC, WB
Recommended Dilution:	WB: 1:200-400; ICC: 1:100-200; IHC: 1:100-200; ELISA: 1:5,000-20,000
Reactivity:	Rabbit IgG
Host:	Mouse
Immunogen:	Rabbit IgG, whole molecule
Isotype:	IgG
Formulation:	0.01M PBS (pH7.4), 0.01% Thimerosal, 50% glycerol
Concentration:	lot specific
Purification:	This antibody was purified from antiserum by immunoaffinity chromatography which removes essentially all mouse serum proteins, except the specific antibody for rabbit IgG.
Conjugation:	Biotin
Storage:	Store at -20°C as received. Avoid repeated freezing and thawing. After first thaw, pipette into single aliquote for future use.
Note:	This biotin conjugated antibody is specific for rabbit IgG and shows no cross-reactivity with rabbit IgA, IgM and other animal species IgG proteins.



[View online »](#)