

Product datasheet for **TA130021**

FITC Conjugated Goat Anti-Rabbit IgG

Product data:

Product Type:	Secondary Antibodies
Product Name:	FITC Conjugated Goat Anti-Rabbit IgG
Applications:	FC, IHC
Recommended Dilution:	FC: 1:100; IF/ICC: 1:32-64; IHC: 1:32-64
Reactivity:	Rabbit IgG
Host:	Goat
Immunogen:	Rabbit IgG, whole molecule
Isotype:	IgG
Formulation:	0.01M PBS (pH7.4), 0.01% Thimerosal, 50% glycerol
Concentration:	lot specific
Purification:	This antibody was purified from antiserum by immunoaffinity chromatography which removes essentially all goat serum proteins, except the specific antibody for rabbit IgG.
Conjugation:	FITC
Storage:	Store at -20°C as received. Avoid repeated freezing and thawing. After first thaw, pipette into single aliquote for future use.
Note:	This FITC conjugated antibody is specific for rabbit IgG and shows no cross-reactivity with rabbit IgA, IgM and other animal species IgG proteins.



[View online »](#)