

Product datasheet for **TA130024**

Peroxidase Conjugated Goat Anti-Rabbit IgG (gamma-chain specific)

Product data:

Product Type:	Secondary Antibodies
Product Name:	Peroxidase Conjugated Goat Anti-Rabbit IgG (gamma-chain specific)
Applications:	ELISA, WB
Recommended Dilution:	WB: 1:3,000-10,000; ELISA: 1:20,000-40,000
Reactivity:	Rabbit IgG
Host:	Goat
Immunogen:	Rabbit IgG, whole molecule
Isotype:	IgG
Formulation:	0.01M PBS (pH7.4), 50% glycerol
Concentration:	lot specific
Purification:	This antibody was purified from antiserum by immunoaffinity chromatography which removes essentially all goat serum proteins, except the specific antibody for rabbit IgG.
Conjugation:	HRP
Storage:	Store at -20°C as received. Avoid repeated freezing and thawing. After first thaw, pipette into single aliquote for future use.
Note:	This peroxidase conjugated antibody is specific for rabbit IgG and shows no cross-reactivity with rabbit IgA, IgM and other animal species IgG proteins.



[View online »](#)