

Product datasheet for **TA130025**

Rabbit Anti-Goat IgG Secondary Antibody

Product data:

Product Type:	Secondary Antibodies
Product Name:	Rabbit Anti-Goat IgG Secondary Antibody
Applications:	ELISA, IHC, WB
Recommended Dilution:	WB, IHC, IF/ICC, ELISA
Reactivity:	Goat IgG
Host:	Rabbit
Immunogen:	Goat IgG, whole molecule
Isotype:	IgG
Formulation:	Lyophilized with 1% sodium acetate
Purification:	This antibody was purified from antiserum by immunoaffinity chromatography which removes essentially all rabbit serum proteins, except the specific antibody for goat IgG.
Conjugation:	Unconjugated
Storage:	Store at -20°C as received. Avoid repeated freezing and thawing. After first thaw, pipette into single aliquote for future use.



[View online »](#)