

Product datasheet for **TA130028**

Cy3 Conjugated Rabbit Anti-Goat IgG

Product data:

Product Type:	Secondary Antibodies
Product Name:	Cy3 Conjugated Rabbit Anti-Goat IgG
Applications:	IHC
Recommended Dilution:	IF/ICC: 1:50-100; IHC: 1:50-100
Reactivity:	Goat IgG
Host:	Rabbit
Immunogen:	Goat IgG, whole molecule
Isotype:	IgG
Formulation:	0.01M PBS (pH7.4), 0.01% Thimerosal, 50% glycerol
Concentration:	lot specific
Purification:	This antibody was purified from antiserum by immunoaffinity chromatography which removes essentially all rabbit serum proteins, except the specific antibody for goat IgG.
Conjugation:	Cy3
Storage:	Store at -20°C as received. Avoid repeated freezing and thawing. After first thaw, pipette into single aliquote for future use.



[View online »](#)