

## Product datasheet for **TA130030**

### **FITC Conjugated Rabbit Anti-Goat IgG**

#### **Product data:**

<b>Product Type:</b>	Secondary Antibodies
<b>Product Name:</b>	FITC Conjugated Rabbit Anti-Goat IgG
<b>Applications:</b>	FC, IHC
<b>Recommended Dilution:</b>	FC: 1:100; IF/ICC: 1:32-64; IHC: 1:32-64
<b>Reactivity:</b>	Goat IgG
<b>Host:</b>	Rabbit
<b>Immunogen:</b>	Goat IgG, whole molecule
<b>Isotype:</b>	IgG
<b>Formulation:</b>	0.01M PBS (pH7.4), 0.01% Thimerosal, 50% glycerol
<b>Concentration:</b>	lot specific
<b>Purification:</b>	This antibody was purified from antiserum by immunoaffinity chromatography which removes essentially all rabbit serum proteins, except the specific antibody for goat IgG.
<b>Conjugation:</b>	FITC
<b>Storage:</b>	Store at -20°C as received. Avoid repeated freezing and thawing. After first thaw, pipette into single aliquote for future use.
<b>Note:</b>	This FITC conjugated antibody is specific for goat IgG and shows no cross-reactivity with goat IgA, IgM and other animal species IgG proteins.



[View online »](#)