

Name: Rabbit Anti-Rat IgG Secondary Antibody
Product Data Sheet

Catalog: TA130033

Target:	Rat IgG
Immunogen:	Rat IgG, whole molecule
Conjugation:	Unconjugated
Concentration:	
Host Species:	Rabbit
IsoType:	IgG
Application:	WB, IHC, IF/ICC, ELISA
Form:	Lyophilized
Buffer:	Lyophilized with 1% sodium acetate
Purification:	This antibody was purified from antiserum by immunoaffinity chromatography which removes essentially all rabbit serum proteins, except the specific antibody for rat IgG.
Storage:	At -20°C for 2 years. Avoid repeated freezing and thawing.

Note:

* More validation images may be available on our website:
http://www.origene.com/secondary_antibody/TA130033.aspx

Related Product:

- TrueORF cDNA clones
- VERIFY Tagged Antigen lysates
- HuSH-29 shRNA
- Western Blot reagents
- Anti-myc/DDK tag antibodies

This product is to be used for laboratory only. Not for diagnostic or therapeutic use.

© 2017 OriGene Technologies, Inc. 9620 Medical Center Dr., Suite 200, Rockville, MD 20850
