

## Product datasheet for **TA130045**

### FITC Conjugated Goat Anti-human IgA

#### Product data:

|                       |                                                                                                                                                                           |
|-----------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Type:         | Secondary Antibodies                                                                                                                                                      |
| Product Name:         | FITC Conjugated Goat Anti-human IgA                                                                                                                                       |
| Applications:         | FC, IHC                                                                                                                                                                   |
| Recommended Dilution: | FC: 1:100; IF/ICC: 1:32-64; IHC: 1:32-64                                                                                                                                  |
| Reactivity:           | Human IgA                                                                                                                                                                 |
| Host:                 | Goat                                                                                                                                                                      |
| Immunogen:            | Human IgA, whole molecule                                                                                                                                                 |
| Isotype:              | IgA                                                                                                                                                                       |
| Formulation:          | 0.01M PBS (pH7.4), 0.01% Thimerosal, 50% glycerol                                                                                                                         |
| Concentration:        | lot specific                                                                                                                                                              |
| Purification:         | This antibody was purified from antiserum by immunoaffinity chromatography which removes essentially all goat serum proteins, except the specific antibody for human IgA. |
| Conjugation:          | FITC                                                                                                                                                                      |
| Storage:              | Store at -20°C as received. Avoid repeated freezing and thawing. After first thaw, pipette into single aliquote for future use.                                           |
| Note:                 | This FITC conjugated antibody is specific for human IgA and shows no cross-reactivity with human IgG, IgM and other animal species IgA proteins.                          |



[View online »](#)