

Product datasheet for **TA130046**

FITC Conjugated Rabbit Anti-human IgG

Product data:

Product Type:	Secondary Antibodies
Product Name:	FITC Conjugated Rabbit Anti-human IgG
Applications:	FC, IHC
Recommended Dilution:	FC: 1:100; IF/ICC: 1:32-64; IHC: 1:32-64
Reactivity:	Human IgG
Host:	Rabbit
Immunogen:	Human IgG, whole molecule
Isotype:	IgG
Formulation:	0.01M PBS (pH7.4), 0.01% Thimerosal, 50% glycerol
Concentration:	lot specific
Purification:	This antibody was purified from antiserum by immunoaffinity chromatography which removes essentially all rabbit serum proteins, except the specific antibody for human IgG.
Conjugation:	FITC
Storage:	Store at -20°C as received. Avoid repeated freezing and thawing. After first thaw, pipette into single aliquote for future use.
Note:	This FITC conjugated antibody is specific for human IgG and shows no cross-reactivity with human IgA, IgM and other animal species IgG proteins.



[View online »](#)