

Product datasheet for **TA130049**

Peroxidase Conjugated Goat Anti-human IgA (alpha-chain specific)

Product data:

Product Type:	Secondary Antibodies
Product Name:	Peroxidase Conjugated Goat Anti-human IgA (alpha-chain specific)
Applications:	ELISA, WB
Recommended Dilution:	WB: 1:3,000-10,000; ELISA: 1:20,000-40,000
Reactivity:	Human IgA
Host:	Goat
Immunogen:	Human IgA, whole molecule
Isotype:	IgA
Formulation:	0.01M PBS (pH7.4), 50% glycerol
Concentration:	lot specific
Purification:	This antibody was purified from antiserum by immunoaffinity chromatography which removes essentially all goat serum proteins, except the specific antibody for human IgA.
Conjugation:	HRP
Storage:	Store at -20°C as received. Avoid repeated freezing and thawing. After first thaw, pipette into single aliquote for future use.
Note:	This peroxidase conjugated antibody is specific for human IgA and shows no cross-reactivity with human IgG, IgM and other animal species IgA proteins.



[View online »](#)