

Product datasheet for **TA130060**

Goat Monoclonal Anti-Rabbit IgG Fc antibody, clone RMG02

Product data:

Product Type:	Secondary Antibodies
Product Name:	Goat Monoclonal Anti-Rabbit IgG Fc antibody, clone RMG02
Applications:	ELISA, WB
Recommended Dilution:	ELISA: 0.01ug/mL – 0.5ug/mL, ICC/IHC: 0.5ug/mL-2ug/mL, WB (non-reduced): 0.1ug/mL-0.5ug/mL
Reactivity:	Rabbit IgG Fc
Host:	Goat
Immunogen:	Rabbit IgG
Isotype:	IgG
Formulation:	PBS with 1% BSA and 0.09% sodium azide
Concentration:	lot specific
Purification:	Protein A affinity purified from an animal origin-free culture supernatant
Conjugation:	Unconjugated
Storage:	Store at -20°C as received. Avoid repeated freezing and thawing. After first thaw, pipette into single aliquote for future use.
Note:	This antibody reacts to the Fc region of Rabbit IgG. No cross reactivity with human IgG, rat IgG, or mouse IgG.



[View online »](#)