

## Product datasheet for **TA130061**

### Rabbit Monoclonal Anti-Human Gamma Heavy Chain antibody, clone RMG02

#### Product data:

Product Type:	Secondary Antibodies
Product Name:	Rabbit Monoclonal Anti-Human Gamma Heavy Chain antibody, clone RMG02
Applications:	ELISA, ICC, IHC
Recommended Dilution:	ELISA: 0.05ug/mL – 0.2ug/mL, ICC, IHC: 0.5ug/mL-2ug/mL
Reactivity:	Human gamma heavy chain
Host:	Rabbit
Immunogen:	Human IgG
Isotype:	IgG
Formulation:	PBS with 1% BSA and 0.09% sodium azide
Concentration:	lot specific
Purification:	Protein A affinity purified from an animal origin-free culture supernatant
Conjugation:	Unconjugated
Storage:	Store at -20°C as received. Avoid repeated freezing and thawing. After first thaw, pipette into single aliquote for future use.
Note:	This antibody reacts to the Fc region of all gamma heavy chains of human immunoglobulins, including $\gamma 1$ , $\gamma 2$ , $\gamma 3$ , and $\gamma 4$ . No cross reactivity with other human heavy chains, mouse IgG, rat IgG or goat IgG



[View online »](#)