

Product datasheet for **TA130062**

Rabbit Monoclonal Anti-Human Kappa Light Chain antibody, clone RM126

Product data:

Product Type:	Secondary Antibodies
Product Name:	Rabbit Monoclonal Anti-Human Kappa Light Chain antibody, clone RM126
Applications:	ELISA
Recommended Dilution:	ELISA: 0.05ug/mL-0.2ug/ml; ICC/IHC: 0.5ug/mL-2ug/mL
Reactivity:	Human kappa light chain
Host:	Rabbit
Immunogen:	Human IgG
Isotype:	IgG
Formulation:	PBS with 1% BSA and 0.09% sodium azide
Concentration:	lot specific
Purification:	Protein A affinity purified from an animal origin-free culture supernatant
Conjugation:	Unconjugated
Storage:	Store at -20°C as received. Avoid repeated freezing and thawing. After first thaw, pipette into single aliquote for future use.
Note:	This antibody reacts to the kappa light chain of human immunoglobulins. No cross reactivity with the lambda light chain, mouse IgG, rat IgG, or goat IgG.



[View online »](#)