

Product datasheet for **TA130063**

Rabbit Monoclonal Anti-Human Lambda Light Chain antibody, clone RM127

Product data:

Product Type:	Secondary Antibodies
Product Name:	Rabbit Monoclonal Anti-Human Lambda Light Chain antibody, clone RM127
Applications:	ELISA
Recommended Dilution:	ELISA: 0.05ug/mL-0.2ug/ml; ICC/IHC: 0.5ug/mL-2ug/mL
Reactivity:	Human lambda light chain
Host:	Rabbit
Immunogen:	Human IgG
Isotype:	IgG
Formulation:	PBS with 1% BSA and 0.09% sodium azide
Concentration:	lot specific
Purification:	Protein A affinity purified from an animal origin-free culture supernatant
Conjugation:	Unconjugated
Storage:	Store at -20°C as received. Avoid repeated freezing and thawing. After first thaw, pipette into single aliquote for future use.
Note:	This antibody reacts to the lambda light chain of human immunoglobulins. No cross reactivity with the kappa light chain, mouse IgG, rat IgG, or goat IgG



[View online »](#)