

## Product datasheet for **TA130067**

### Rabbit Monoclonal Anti-Human IgA1 antibody, clone RM124

#### Product data:

Product Type:	Secondary Antibodies
Product Name:	Rabbit Monoclonal Anti-Human IgA1 antibody, clone RM124
Applications:	ELISA, ICC, IHC
Recommended Dilution:	ELISA: 50ng/well – 200ng/well (for Capture); 0.05ug/mL – 0.2ug/mL (for Detection); IHC, ICC: 0.5ug/mL – 2ug/mL
Reactivity:	Human IgA1
Host:	Rabbit
Immunogen:	Human IgA
Isotype:	IgG
Formulation:	PBS with 1% BSA and 0.09% sodium azide
Concentration:	lot specific
Purification:	Protein A affinity purified from an animal origin-free culture supernatant
Conjugation:	Unconjugated
Storage:	Store at -20°C as received. Avoid repeated freezing and thawing. After first thaw, pipette into single aliquote for future use.
Note:	This antibody reacts to human IgA1, and very slightly cross reacts with IgA2. No cross reactivity with human IgG, IgM, IgD, or IgE.



[View online »](#)