

Product datasheet for **TA130069**

Biotin Conjugated Rabbit Monoclonal Anti-Human Ig Light Chains antibody 3, clone RM129

Product data:

Product Type:	Secondary Antibodies
Product Name:	Biotin Conjugated Rabbit Monoclonal Anti-Human Ig Light Chains antibody 3, clone RM129
Applications:	ELISA, ICC, IHC
Recommended Dilution:	ELISA: 0.02ug/mL – 0.25ug/mL; IHC, ICC: 0.5ug/mL – 2ug/mL
Reactivity:	Human Ig light chain
Host:	Rabbit
Immunogen:	Human IgG
Isotype:	IgG
Formulation:	PBS with 1% BSA and 0.09% sodium azide
Concentration:	lot specific
Purification:	Protein A affinity purified from an animal origin-free culture supernatant
Conjugation:	Biotin
Storage:	Store at -20°C as received. Avoid repeated freezing and thawing. After first thaw, pipette into single aliquote for future use.
Note:	This antibody reacts to human IgA, including both IgA1 and IgA2. No cross reactivity with human IgG, IgM, IgD, or IgE.



[View online »](#)