

## Product datasheet for **TA130070**

### Goat Monoclonal Anti-Mouse IgG Fab Antibody, Clone RMG05

#### Product data:

Product Type:	Secondary Antibodies
Product Name:	Goat Monoclonal Anti-Mouse IgG Fab Antibody, Clone RMG05
Applications:	ELISA, ICC, IHC, WB
Recommended Dilution:	ELISA: 0.05ug/mL - 1ug/mL; ICC, IHC: 0.5ug/mL - 2ug/mL; WB: 0.1ug/mL - 0.5ug/mL
Reactivity:	Mouse IgG
Host:	Goat
Immunogen:	Mouse IgG
Isotype:	IgG
Formulation:	PBS with 1% BSA and 0.09% sodium azide
Concentration:	lot specific
Purification:	Protein G affinity purified from an animal origin-free culture supernatant
Conjugation:	Unconjugated
Storage:	Store at -20°C as received. Avoid repeated freezing and thawing. After first thaw, pipette into single aliquote for future use.
Note:	This antibody reacts to the Fab region of Mouse IgG, including IgG1, IgG2a, IgG2b, and IgG3. No cross reactivity with IgM, IgA, IgE, human IgG, rat IgG, or rabbit IgG.



[View online »](#)