

Product datasheet for **TA130071**

Goat Monoclonal Anti-Mouse IgG Fc Antibody, Clone RMG06

Product data:

Product Type:	Secondary Antibodies
Product Name:	Goat Monoclonal Anti-Mouse IgG Fc Antibody, Clone RMG06
Applications:	ELISA, ICC, IHC, WB
Recommended Dilution:	ELISA: 0.05ug/mL - 1ug/mL; ICC, IHC: 0.5ug/mL - 2ug/mL; WB: 0.1ug/mL - 0.5ug/mL
Reactivity:	Mouse IgG
Host:	Goat
Immunogen:	Mouse IgG
Isotype:	IgG
Formulation:	PBS with 1% BSA and 0.09% sodium azide
Concentration:	lot specific
Purification:	Protein G affinity purified from an animal origin-free culture supernatant
Conjugation:	Unconjugated
Storage:	Store at -20°C as received. Avoid repeated freezing and thawing. After first thaw, pipette into single aliquote for future use.
Note:	This antibody reacts to the Fc region of Mouse IgG, including IgG1, IgG2a, and IgG2b. It shows very slight reaction with mouse IgG3. No cross reactivity with IgM, IgA, IgE, human IgG, rat IgG, or rabbit IgG.



[View online »](#)