

Product datasheet for **TA130076**

Rabbit Monoclonal Anti-Human IgD Antibody, Clone RM123

Product data:

Product Type:	Secondary Antibodies
Product Name:	Rabbit Monoclonal Anti-Human IgD Antibody, Clone RM123
Applications:	ELISA, ICC, IHC
Recommended Dilution:	ELISA: 25ng/well - 200ng/well for Capture;; 0.01 ug/mL - 0.1ug/mL for Detection;; , ICC: 0.5ug/mL-2ug/mL, IHC: 0.5ug/mL-2ug/mL
Reactivity:	Human IgD
Host:	Rabbit
Immunogen:	Human IgD
Isotype:	IgG
Formulation:	PBS with 1% BSA and 0.09% sodium azide
Concentration:	lot specific
Purification:	Protein A affinity purified from an animal origin-free culture supernatant
Conjugation:	Unconjugated
Storage:	Store at -20°C as received. Avoid repeated freezing and thawing. After first thaw, pipette into single aliquote for future use.
Note:	RM123 reacts to human IgD. No cross reactivity with human IgG, IgM, IgE, or IgA.



[View online »](#)