

Product datasheet for **TA150030**

Clone OTI4C5, Anti-DDK (FLAG) Monoclonal antibody, HRP


Product data:

Product Type:	Tag Antibodies
Product Name:	Clone OTI4C5, Anti-DDK (FLAG) Monoclonal antibody, HRP
Clone Name:	OTI4C5 (formerly 4C5)
Applications:	IP, WB
Recommended Dilution:	WB: 1:1000-1:5000, IP: 1:50
Reactivity:	DDK (FLAG tag): DYKDDDDK
Host:	Mouse
Clonality:	Monoclonal
Immunogen:	Anti-DDK monoclonal antibody is produced by immunizing mice with a synthetic peptide coupled to KLH
Isotype:	IgG2a
Formulation:	PBS (pH7.4) containing 1% BSA and 40% glycerol
Concentration:	0.5mg/ml
Conjugation:	HRP
Storage:	Store at -20°C. For maximum recovery of product, centrifuge the original vial after thawing and prior to opening the cap. Aliquot to avoid repeated freezing and thawing.
Stability:	Stable for 1 year at -20°C from date of shipment.
Synonyms:	ANTI DDK, ANTI-DDK, ANTIDDK, FLAG TAG, ANTI FLAG, ANTI-FLAG, ANTI FLAG ANTIBODY, FLAG ANTIBODY, M2 ANTIBODY, FLAG M2, ANTIBODY FOR FLAG, 4C5, ANTI DDK 4C5, DDK 4C5, ANTI-DDK 4C5, DYKDDDDK TAG, D-TAG ANTIBODY, ANTI D-TAG, ANTI D TAG, ANTI-D TAG, ANTI-MDDK, ANTI-MYC DDK, ANTI-MYC/DDK, MYC/DDK ANTIBODY, MYC DDK ANTIBODY, ANTI MYC DDK, ANTI DDK MYC, ANTI-DDK/MYC, ANTI-DDK MYC, FLAG AB, FLAG MAB, DDK ANTIBODY, ANTI FLAG HRP, ANTI DDK HRP, ANTI-FLAG HRP, ANTI-DDK HRP
Note:	This antibody doesn't recognize N- terminus DDK tag. Please use TA180144 for the detection of DDK tag at any location in Western blotting.



[View online »](#)

Product images:



Myc/DDK tagged VARS2 overexpression lysate and pull-down by various batches of anti-DDK beads were detected with HRP-conjugated 4C5 anti-DDK antibody (TA150030).