

Product datasheet for **TA150033**

Anti-GST rabbit polyclonal antibody


Product data:

Product Type:	Tag Antibodies
Product Name:	Anti-GST rabbit polyclonal antibody
Applications:	IP, WB
Recommended Dilution:	WB: 1:3000, IP: 1:3000
Reactivity:	GST tag
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	Rabbits were immunized with recombinant Glutathione-S-transferase (GST) from E.coli . After multiple immunizations in Freund's adjuvant serum was collected. Antibodies were immunoaffinity purified using GST immobilized on a solid support.
Isotype:	IgG
Formulation:	PBS (pH7.0) containing 1% BSA, 10% glycerol
Concentration:	lot specific
Conjugation:	Unconjugated
Storage:	Store at -20°C. For maximum recovery of product, centrifuge the original vial after thawing and prior to opening the cap. Aliquot to avoid repeated freezing and thawing.
Stability:	Stable for 1 year at -20°C from date of shipment.
Synonyms:	ANTI-GST; ANTIBODY FOR GST; ANTI GST; GST AB; GST ANTIBODY




[View online »](#)

Product images:



Sample: A: 20ug of E. coli only; B: 20ug of Transformed E. coli; C: 20ug of Transformed E. coli and IPTG induced. 12% SDS PAGE, TA150033 diluted at 1:3000.



GST-tagged protein (Mw: 45 kDa) was immunoprecipitated from E. coli lysate using 8 ug of anti-GST antibody (TA150033). The precipitated protein was detected by TA150033 diluted at 1:3000.