

Product datasheet for **TA150043**


Mouse monoclonal tGFP antibody (clone OTI2H8), HRP conjugated

Product data:

Product Type:	Tag Antibodies
Product Name:	Mouse monoclonal tGFP antibody (clone OTI2H8), HRP conjugated
Clone Name:	OTI2H8 (formerly 2H8)
Applications:	WB
Recommended Dilution:	WB: 1:1000
Target:	tGFP
Host:	Mouse
Clonality:	Monoclonal
Immunogen:	Anti-tGFP monoclonal antibody was produced by immunizing mice with purified tGFP protein expressed in HEK293T cell.
Isotype:	IgG2b
Formulation:	PBS (pH7.4) containing 0.1% BSA, 50% Glycerol with final concentration 500ug/ml
Concentration:	0.5mg/ml
Conjugation:	HRP
Storage:	Store at -20°C. For maximum recovery of product, centrifuge the original vial after thawing and prior to opening the cap. Aliquot to avoid repeated freezing and thawing.
Stability:	Stable for 1 year at -20°C from date of shipment.
Predicted Protein Size:	26 kDa



[View online »](#)

Product images:

HEK293T cells were transfected with pCMV6-AC-tGFPvector (Cat# [PS100010], Lane 1) or pCMV6-Entry-MCR2 (Myc-DDK-tagged ORF clone of human MCR2 gene, Lane 2) plasmids for 48 hrs. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-tGFP-HRP (1:1000).