

## Product datasheet for TA150046

### Loading Control Antibody Sample Kit

#### Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB 1:500-1:2000
Reactivity:	Dog, Human, Monkey, Mouse, Rat
Host:	Mouse
Clonality:	Monoclonal
Formulation:	TBS (pH7.4) containing 0.1% BSA, 0.05% NaN <sub>3</sub> and 50% Glycerol
Concentration:	1 mg/ml
Purification:	Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)
Conjugation:	Unconjugated
Storage:	Store at -20°C as received
Stability:	Stable for 12 months from date of receipt.



[View online »](#)