

Product datasheet for **TA150052**

Mouse monoclonal eGFP antibody, clone OT15A2


Product data:

Product Type:	Tag Antibodies
Product Name:	Mouse monoclonal eGFP antibody, clone OT15A2
Clone Name:	OT15A2 (formerly 5A2)
Applications:	WB
Recommended Dilution:	WB 1:500-1:2000
Reactivity:	eGFP
Host:	Mouse
Clonality:	Monoclonal
Immunogen:	A full length eGFP tagged recombinant protein expressed in HEK293T cell
Isotype:	IgG1
Formulation:	PBS (pH 7.3) containing 1%BSA, 50% glycerol and 0.02% sodium azide
Concentration:	0.7mg/ml
Conjugation:	Unconjugated
Storage:	Store at -20°C.
Stability:	Stable for 1 year at -20°C from date of shipment.
Predicted Protein Size:	26.9 kDa



[View online »](#)

Product images:



HEK293T cells were transfected with the pCMV6-ENTRY control (left lane) or pCMV6-ENTRY eGFP (right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immuno-blotted with 5A2 anti-eGFP monoclonal antibody (TA150052) (1:2000).




Figure from citation: Western Blot of eGFP-fused protein level by using anti-eGFP antibody in HeLa cells. Dilution: 1:1000 [View Citation](#)