

## Product datasheet for **TA150061**

### Mouse monoclonal turboRFP antibody, clone OTI2D9


#### Product data:

Product Type:	Tag Antibodies
Product Name:	Mouse monoclonal turboRFP antibody, clone OTI2D9
Clone Name:	OTI2D9 (formerly 2D9)
Applications:	WB
Recommended Dilution:	WB: 1:2000
Reactivity:	tRFP
Host:	Mouse
Clonality:	Monoclonal
Immunogen:	recombinant tRFP produced in E.coil
Isotype:	IgG1
Formulation:	PBS (pH 7.3) containing 1%BSA, 50% glycerol and 0.02% sodium azide
Concentration:	1mg/ml
Conjugation:	Unconjugated
Storage:	Store at -20°C.
Stability:	Stable for 1 year at -20°C from date of shipment.
Predicted Protein Size:	26 kDa



[View online »](#)

## Product images:



Anti-tRFP immunoblot analysis of HEK293T over-expression lysates containing recombinant tRFP, TurboFP602 TurboFP635 fluorescence proteins at 1:2000 dilution.