


## Product datasheet for TA150089

### Anti-Tdtomato chicken polyclonal antibody

#### Product data:

Product Type:	Tag Antibodies
Product Name:	Anti-Tdtomato chicken polyclonal antibody
Applications:	WB
Recommended Dilution:	WB 1:500-1:2000
Reactivity:	Tdtomato
Host:	Chicken
Clonality:	Polyclonal
Immunogen:	Purified Tdtomato protein expressed in E. coli.
Isotype:	IgY
Formulation:	1XPBS, pH 7.2 with 0.05% Sodium azide
Concentration:	1.0 mg/ml. Crude IgY preparation from egg yolks using a precipitation method (Polyethylene Glycol).
Conjugation:	Unconjugated
Storage:	Store at -20°C. For maximum recovery of product, centrifuge the original vial after thawing and prior to opening the cap. Aliquot to avoid repeated freezing and thawing.
Stability:	Stable for 1 year at -20°C from date of shipment.

#### Product images:



HEK293T cells were transfected with the pCMV6-ENTRY tdTomato cDNA. 2ug cell lysate was separated by SDS-PAGE and immunoblotted with chicken anti-TdTomato polyclonal antibody (TA150089, 1:2000)



[View online »](#)