

Product datasheet for **TA150094M**

Clone OTI7E11, Anti-V5 tag monoclonal antibody


Product data:

Product Type:	Tag Antibodies
Product Name:	Clone OTI7E11, Anti-V5 tag monoclonal antibody
Clone Name:	OTI7E11 (formerly 7E11)
Applications:	WB
Recommended Dilution:	WB: 1:2000
Reactivity:	V5 tag
Host:	Mouse
Clonality:	Monoclonal
Immunogen:	Anti-V5 monoclonal antibody is produced by immunizing mice with a synthetic peptide coupled to KLH
Isotype:	IgG1
Formulation:	0.92 mg/ml in PBS (pH 7.3) containing 1%BSA, 50% glycerol and 0.02% sodium azide
Concentration:	1 mg/ml
Conjugation:	Unconjugated
Storage:	Store at -20°C.
Stability:	Stable for 1 year at -20°C from date of shipment.
Synonyms:	ANTI-V5; ANTI-V5 TAG; ANTI V5 TAG; ANTI V5-TAG; V5 ANTIBODY; V5 TAG; V5 TAG ANTIBODY



[View online »](#)

Product images:



Anti-V5 immunoblot analysis of six different HEK293T over-expression lysates containing different V5 tagged recombinant proteins at 1:2000 dilution.