


Product datasheet for TA150122

Rabbit Polyclonal mGFP Antibody

Product data:

Product Type:	Tag Antibodies
Product Name:	Rabbit Polyclonal mGFP Antibody
Applications:	WB
Recommended Dilution:	WB: 1:5000
Target:	mGFP, eGFP, mYFP
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	mGFP
Isotype:	IgG
Formulation:	Preservative: 0.05% Sodium azide
Concentration:	This antibody is provided as whole antiserum. The exact antibody concentration is hard to determine.
Conjugation:	Unconjugated
Storage:	Store at -20°C. For maximum recovery of product, centrifuge the original vial after thawing and prior to opening the cap. Aliquot to avoid repeated freezing and thawing.
Stability:	Stable for 1 year at -20°C from date of shipment.

Product images:



HEK293T cells were transfected with the pCMV6-ENTRY mGFP cDNA and the cell lysates were collected after 48 hours. 0µg, 0.1µg, 0.2µg, 0.4µg, 0.8µg, 1.6µg of lysates (lane 1 to lane 6 respectively) were separated by SDS-PAGE and immunoblotted with Rabbit mGFP polyclonal antibody (TA150122) at 1:5000.



[View online »](#)