


## Product datasheet for TA150124

### Chicken Polyclonal mGFP Antibody

#### Product data:

Product Type:	Tag Antibodies
Product Name:	Chicken Polyclonal mGFP Antibody
Applications:	WB
Recommended Dilution:	WB: 1:3000
Reactivity:	mGFP, eGFP, mYFP
Host:	Chicken
Clonality:	Polyclonal
Immunogen:	mGFP
Isotype:	IgY
Formulation:	PBS, pH 7.2 with 0.05% Sodium azide
Concentration:	1 mg/ml
Conjugation:	Unconjugated
Storage:	Store at -20°C. For maximum recovery of product, centrifuge the original vial after thawing and prior to opening the cap. Aliquot to avoid repeated freezing and thawing.
Stability:	Stable for 1 year at -20°C from date of shipment.

#### Product images:



HEK293T cells were transfected with the pCMV6-ENTRY mGFP cDNA and the cell lysates were collected after 48 hours. 0ug, 0.5ug, 1ug, 2ug, 4ug, 8ug of lysates (lane 1 to lane 6 respectively) were separated by SDS-PAGE and immunoblotted with Chicken mGFP polyclonal antibody (TA150124) at 1:3000.



[View online »](#)