

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Product datasheet for TA160051

SARS-CoV-2/SARS-CoV Spike Protein Rabbit Polyclonal Antibody


Product data:

Product Type:	Primary Antibodies
Applications:	ELISA, WB
Recommended Dilution:	ELISA: the antibody will detect 2 ng of free peptide at 1 μg/ml. WB: 4 μg/ml. SARS-CoV-2/SARS-CoV Spike antibody can be used for the detection of SARS-CoV-2/SARS-CoV Spike protein in ELISA and in Western Blot.
Reactivity:	SARS-CoV, SARS-CoV-2
Host:	Rabbit
lsotype:	IgG
Clonality:	Polyclonal
Immunogen:	The antibody was raised against a peptide corresponding to 16 amino acids between aa 410- 460 near the center of SARS-CoV Spike glycoprotein.
Formulation:	PBS containing 0.02% sodium azide.
Concentration:	1 mg/ml
Purification:	Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Database Link:	<u>P59594</u>
Background:	SARS Spike Antibody: A novel coronavirus has recently been identified as the causative agent of SARS (Severe Acute Respiratory Syndrome). Coronaviruses are a major cause of upper respiratory diseases in humans. The genomes of these viruses are positive-stranded RNA approximately 27-31kb in length. SARS infection can be mediated by the binding of the viral spike protein, a glycosylated 139 kDa protein and the major surface antigen of the virus, to the angiotensin-converting enzyme 2 (ACE2) on target cells. This binding can be blocked by a soluble form of ACE2.




This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2023 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US

Product images:



Western Blot Validation with Human Recombinant Protein Loading: 50 ng per lane of human SARS-CoV-2 Spike S1 recombinant protein. Antibodies: SARS-CoV-2/SARS-CoV Spike 3221 (4 µg/mL), 1h incubation at RT in 5% NFDM/TBST. Secondary: Goat anti-rabbit IgG HRP conjugate at 1:10000 dilution.



ELISA Test

Antibodies: SARS-CoV Spike Antibody, TA160051 (1 µg/mL). A direct ELISA was performed using immunogen as coating antigen and the anti-SARS-CoV Spike antibody as the capture antibody. Secondary: Goat anti-rabbit IgG HRP conjugate at 1:20000 dilution. Detection range is from 2 ng/mL to 1000ng/mL.

This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2023 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US